

Procedure for Artemis Dynamic Air Mattress delivery and collection- Covid 19 Wards

Delivery




- Mattresses will be delivered in a CLEAR bag containing the red collection bag and decontamination form- please place these items on the bed frame under the mattress as they will be needed for bagging procedure later.
- Mattresses will be delivered to the ward entrance by DDV technicians (or porters out of hours) and **must be collected by ward staff**



Bagging up at end of use for Collection

- All mattresses must be double bagged & sealed with either a bag tie or tape
- A completed decontamination form must be placed either in the second red bag so it is visible or taped securely to the outside of the bag
- All red bagged mattresses must be placed outside the ward doors for collection by DDV technicians and must be separated from clean mattresses
- **FAILURE TO COMPLY WITH THIS PROCEDURE MEANS THAT YOUR MATTRESSES WILL NOT BE COLLECTED**



			Sidhill Business Park, Holmfield, Halifax, West Yorkshire, H10 2TN Telephone: 0800 0370234 Email: service@driveandvilbiss.co.uk Web: www.sidhill.co.uk
DECONTAMINATION CERTIFICATE			
<small>This is to provide confirmation that this equipment has been cleaned and decontaminated in accordance with Department of Health NHS Guidelines - HTM 01-04</small>			
DATE:	_____		
CUSTOMER – TRUST/WARD/LOCATION:	_____		
TYPE OF PRODUCT (MATTRESS/BED):	_____		
MODEL/MANUFACTURER:	_____		
ID NO'S - SERIAL/RENTAL:	_____		
SIGN & PRINT NAME:	_____		
Rev 02			12 th September 2017

* DDV= Drive DeVilbiss (previously Sidhill)